


Beltone **Legend™**

Your perfect match

Personalized hearing



Made for
iPod iPhone iPad

 **Beltone**



Taking **personalized hearing** to the next level

Every person is unique, and so is their hearing. So a hearing aid doesn't just need to match your hearing loss, it needs to match your lifestyle and individual tastes. You know what fits - ergonomically and aesthetically. You pay attention to certain sounds, and ignore others.

Beltone Legend hearing aids adapt to your unique situation and individual needs better than any other hearing aid on the market. Superb sound quality and speech understanding are the main focus, thanks to Personal Sound ID™ and Crosslink Directionality, but

other features like Smart Gain Pro and Feedback Eraser add adaptability, and intelligence. Made for iPhone® integration lets you stream sounds wirelessly to your hearing aids from your iPhone, iPad® or iPod touch®, while the HearPlus app allows you to control your hearing aids directly from your iPhone, iPad and iPod touch and the most popular Android™ phone model*. There is no need to compromise when you choose Beltone Legend.

Beltone Legend really is your perfect match!

Choose your own sound

We're always listening – for voices, ringtones, signals warning of danger. We're selective. Our brains compare sounds and build a three-dimensional sound picture. With most other hearing aids all the subtlety disappears, meaning that you lose your ability to distinguish the specific sound environment.

Restores more natural hearing

Beltone Legend helps restore natural hearing with groundbreaking features: **Personal Sound ID™** and **CrossLink Directionality with Personal Sound ID**. Personal Sound ID mimics the natural sound level differences between your ears, helping you identify where sounds are coming from.

CrossLink Directionality enhances sound quality, monitoring the sound environment and automatically switching to the most natural speech focus settings. CrossLink Directionality with Personal Sound ID delivers sound the way nature intended. It lets you instinctively choose what sounds you want to focus on, while other hearing aids try to artificially make those choices for you.

With these advanced features, focusing on a friend's words becomes easy again, or switching focus to an alarm or cry for help.



MY PERFECT MATCH

Name: Frank

Profession: Winemaker and organizer of wine tastings

Passion: Wine and hiking in the mountains

Hearing loss: Moderate

Frank's choice: Beltone Legend 63DW for discreetness, connectivity to iPhone and maximum speech understanding.

Adapts to every lifestyle

The advanced features in Beltone Legend hearing aids adapt to all kinds of environments, situations and needs – from a quiet chat in the corner to a windy boat ride.

Smart Gain Pro™ adjusts volume, adapting automatically to changes in the environment. With low background noise, little volume is needed. With a lot of ambient noise, volume should be higher.

Sound Cleaner Pro™ works with Smart Gain Pro to minimize background noise so you can move into noisy environments and still hear the whole conversation.

Feedback Eraser™ with WhistleStop™ senses feedback, cancelling it out with no quality loss. It eliminates unwanted noises and whistling when you use the phone, or hug your relatives.

Ear to Ear Synchronization automatically applies program and volume changes from one hearing aid to the other, so you never have to change settings on both instruments.

Sound Shifter™ makes it possible to hear higher-pitched sounds, so you don't miss out on the full spectrum.



Made to fit you

Something for everyone

Beltone Legend stands out as the most flexible hearing aid on the market. They are small and virtually invisible and built for everyday life, from indoors to extreme sports and tough weather. Thanks to numerous details, Beltone Legend is at home anywhere. HPF⁶⁰ NanoBlock coating blocks moisture and dirt. The Wind Noise Reduction feature makes being outdoors in windy weather easy for you, filtering out any uncomfortable noise.

Beltone Legend hearing aids are available in a variety of discreet models and colors to suit any taste.



Invisible-in-Canal (IIC)



Receiver-in-Ear (RIE)



IIC, CIC, ITC, ITE, MIH-S, MIH, 63 RIE, 64 RIE, 76 BTE, 86 BTE



MY PERFECT MATCH

Name: Linda

Profession: Architect

Passion: Loves to be active and travel the world

Hearing loss: Mild

Linda's choice: Beltone Legend IIC for invisibility and custom design.

Control through popular mobile devices

Personal control and direct streaming

Thanks to Beltone Legend's Made for iPhone integration, you can stream sound directly from an iPhone, iPad or iPod touch to your hearing aids without having to wear any sort of adapter around your neck. Stream voices from phone calls, use FaceTime®, listen to music and much more.

Free HearPlus™ app

Download the HearPlus™ app and start controlling your hearing aids in a simple and intuitive way directly from your iPhone, iPad or iPod touch – or the most popular Android phones. Adjust volume or program settings. Use the My Places function for quick custom adjustment to favorite locations or use the 'Finder' function to locate misplaced hearing aids. You can even adjust background noise levels, speech focus or the amount of wind noise reduction directly with the app.



The Beltone HearPlus app



MY PERFECT MATCH

Name: John

Profession: Retired Construction Manager

Passion: Sailing and spending time with his grandchildren

Hearing loss: Severe

John's choice: Beltone Legend 64DW for extended battery life, connectivity to iPhone and wind noise reduction.

Beltone. Helping the world hear better

Beltone was founded on the act of helping a friend to enjoy life more. Since 1940, we have provided knowledge, tools, service and training to the professionals servicing the hearing impaired. Beltone makes quality hearing care available to more people. To do so, we develop technically optimal hearing solutions without ever losing focus on the individual's needs. It is our belief that no individual should be denied the possibility of an improved life with better hearing.

*Beltone Legend is compatible with iPhone 6, iPhone 6 Plus, iPhone 5s, iPhone 5c, iPhone 5, iPad Air 2, iPad Air, iPad (4th generation), iPad mini 3, iPad mini 2, iPad mini and iPod touch (5th generation) using iOS 7.X or later. Apple, the Apple logo, iPhone, iPad, iPod touch and FaceTime are trademarks of Apple Inc., registered in the U.S. and other countries. App Store is a service mark of Apple Inc.

The Android version of the Beltone HearPlus app runs on selected Android devices. See www.beltone-hearing.com for a list of Android phones that work with the app. Android and Google Play are trademarks of Google Inc.



Worldwide Headquarters

Beltone A/S

Lautrupbjerg 7

DK-2750 Ballerup

Denmark

Tel.: + 45 45 75 11 11

Fax: + 45 45 75 11 19

Web: www.beltone-hearing.com

