






# Solus Max Cable Guide

	<b>Grey Cable</b>	<b>Grey Cable</b>	<b>Battery Pill Cables</b>	<b>Grey Cable</b>
	 <p>Cable # SA04305 Red Cable# SA04304 Blue Black Boot # 19902500 Taupe Boot # 17598200</p>	 <p>Cable # SA04305 Red Cable # SA04304 Blue Flex Strip # 19533900</p>	 <p>Part # 18523700 Part # 18523800</p>	 <p>Kit Part # 19215000</p>
<b>Beltone Trust</b>				
BTE/RIE: 86, 76 & 64	● Black			
BTE/RIE: 66 & 63 Custom: ITE, ITC, MIH, CIC, MIH-S		●		
IIC			●	
<b>Beltone Legend</b>				
BTE/RIE: 86, 76 & 64	● Black			
RIE: 63 Custom: ITE, ITC, MIH, CIC, MIH-S		●		
IIC			●	
<b>Beltone Boost Max &amp; Boost Plus</b>				
BTE: 86	● Black			
Super Power BTE: 95				●
<b>Beltone Ally</b>				
BTE/RIE: 86, 76, 66 & 64	● Black			
All Custom Products		●		
<b>Beltone Origin &amp; Bold</b>				
Origin BTE: 85, 75, 65 Bold: 95	● Taupe			
Origin 3 RIE: 64	● Black			
Origin: All Custom Products		●		